**ASPEL S.A.**

32-080 ZABIERZÓW, os. Sienkiewicza 33  
tel. +48 12 285-22-22, fax +48 12 285-30-30  
[www.aspel.com.pl](http://www.aspel.com.pl)

Elektrokardiograf

**AsCARD B5 ECO Mr GREEN**

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WYDANIE IV

---

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**Operation Manual  
Electrocardiograph  
AsCARD B5 ECO Mr GREEN**

ZABIERZÓW, POLAND  
October, 2004

EDITION IV

## **CONGRATULATIONS**

- We are pleased that you have chosen the AsCARD B5 ECO Mr GREEN electrocardiograph;
  - Choosing the services offered by ASPEL you have proved that you do not approve compromise and you want the highest quality.
  - ASPEL has the pleasure to present the AsCARD B5 ECO Mr GREEN electrocardiograph which is a result of years of surveys and experiences gathered through close contact with our customers. You have chosen high durability and efficiency as these are the features of AsCARD B5 ECO Mr GREEN.
  - Read this Manual carefully as it includes guidelines concerning the safe installation, use and maintenance and several practical instructions how to optimise this apparatus application.
  - Store this Manual for reference.
-

## **1. APPLICATION**

The electrocardiograph AsCARD B5 ECO Mr GREEN is a high quality one and three-channel ECG unit which enables to make electrocardiogram in full 12 leads. The ECG unit maintenance is very simple and the modern keyboard makes the operation much easier. The record is made in manual or automatic mode. The unit is supplied from the internal accumulator. The accumulator capacity makes possible to print about 150 examination in auto mode. Low magnitude and small size are the additional features which enable easy transport. The electrocardiograph is designed for carrying out ECG examinations at the patients home and in all health-care institutions.


1. The person operating the electrocardiograph ought to read the Operation Manual carefully before starting the unit operation.
2. The Operation Manual will help the user to operate and maintain electrocardiograph correctly.
3. The efficient functioning and fulfilling the tasks the user reaches after obeying the Operation Manual descriptions.
4. It is forbidden to use the unit with a damaged charger that do not respond to the one given by manufacturer in technical specification.
5. Every repair ought to be made by Authorised Services.
6. The manufacturer is not responsible for any damages resulting from not obeying the Operation Manual.
7. If the electrocardiograph is used in combination with other equipment assure that:
  - all connections are in compliance with international IEC safety standards.
  - all connected units are in compliance with IEC standardsIf the above standards are not obeyed it can cause the safety hazard for which the manufacturer is not responsible.
8. Electrocardiograph used together with a cardiostimulator or other electric stimulator do not cause any safety hazard of the patient or the user.

## 1. 1. Description of unit.

The device is made up of modern microprocessor technology. The unit has a thermal printer, equipped with a head of high resolution to print the ECG leads. The cover made out of ABS resins and the membrane keyboard facilitates cleaning.

## 1. 2. Picture of unit.

1. - paper box
2. - charger socket
3. - charging indicator
4. - LCD display
5. - keyboard
6. - patient cable socket

## 1. 3. Basic Technical Parameter.

Dimensions: (DxSxW) 255x195x66 mm;  
Weight: <2,0 kg  
Power supply: internal accumulator 12 V/0,7 Ah (accumulator with internal charger is build in the unit). The accumulator is changeable by Authorised Service only. Accumulator capacity enables to perform 100 automatic examinations. Remember that the accumulator is self-discharged, so after long period of non-using the unit, the accumulator has to be charged. The capacity of accumulator is decreases during the period of operation, so the number of examinations may lessen. Charging the accumulator is performed with connecting it to main

	supply 115/230 V, 50/60 Hz.
Working temperature:	+10° C ÷ +40° C
Relative humidity:	25% ÷ 95% (no condensation)
Atmospheric pressure:	700 - 1060 hPa
ECG - signals:	12 standard leads I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6;
Sensitivity:	2,5/5/10/20 mm/mV ± 5%;
Speed:	5/25/50 mm/s ±5%;
Paper:	thermoactive, dustfree, 58 mm;
Frequency range:	0,05 ÷ 150 Hz;
Digital filter:	50 Hz i 35 Hz; antidrift filter;
LCD display:	2 x 24 marks;
Using safety:	- protection degree CF (EN 60601-1), - protection class: II ECG input protected against defibrilating
CMRR	>80dB
Fuses:	
- internal supply	<b>T 3,15 A/L/250 V</b> <i>changeable by Authorised Service only</i>
- mains supply	<b>T 500 mA/L/250 V</b> for 230V (T 1,0 A for 115V) <i>changeable by the user</i>
Functional features:	printout in auto mode 5 sec, signalisation of accumulator voltage reduction, memory of last automatic examination, INOP signalisation, signalisation of accumulator charging
Storage and transport:	
Temperature:	+5 to + 40 °C
Relative humidity:	<80% (non-condesive)
Atmospheric pressure:	700 - 1060 hPa
NOTE !	
The unit is not water resistant. Avoid getting damp and wet !	

## 1. 4. Manufacturer.

**MANUFACTURER:**  
**ASPEL S. A.**  
 os. H. Sienkiewicza 33  
 32-080 Zabierzów / POLAND  
 tel. + 48 12 / 285 22 22, fax + 48 12 / 285 30 30

## 1. 5. Marking of unit.

The marking of unit consists of:

Factory number	Year of production	Option							
<table border="1" style="width: 100%; height: 30px;"> <tr> <td style="width: 33%;"></td> <td style="width: 33%;"></td> <td style="width: 33%;"></td> </tr> </table>				<table border="1" style="width: 100%; height: 30px;"> <tr> <td style="width: 50%;"></td> <td style="width: 50%;"></td> </tr> </table>			<table border="1" style="width: 100%; height: 30px;"> <tr> <td style="width: 50%;"></td> <td style="width: 50%;"></td> </tr> </table>		

## 2. WORKING CONDITIONS

Working conditions of electrocardiograph AsCARD B5 ECO Mr GREEN:

- temperature +10 to + 40 °C
- relative humidity from 25% to 95%
- the air should not be polluted with corrosion ingredients.

## 3. FEATURES

The electrocardiograph AsCARD B5 ECO Mr GREEN allows the recording of 12 standard ECG leads. The following modes are possible :

- 1-channel manual record -  
 this mode allows the recording of selected ECG lead. It is possible to set required record speed, sensitivity, filters on/of. It is possible to print one of 12 ECG leads I, II, III; aVR; aVL; aVF; V1; V2; V3; V4; V5; V6;
- 3-channel manual record -  
 this mode allows the recording of selected ECG lead. It is possible to set required record speed, sensitivity, filters on/off. It is possible to print one of four ECG leads I-II-III, aVR-aVL-aVF, V1-V2-V3, V4-V5-V6.
- automatic record -  
 this mode allows the recording of 5 sec of ECG signal from all 12 leads and its printout. The examination is automatically registered in the unit memory. It enables to print the copy of auto examination when requested. Record parameters are the same

as in manual record. The examination is memorised in internal unit memory.

- copy of automatic record -


this mode allows printout of 12 leads from the last automatic examination. Record parameters are set in the the same way as in manual record, but the filters can't be switched on/off. The examination is memorised in internal unit memory.

### 3. 1. How to prepare AsCARD B5 ECO Mr GREEN

The operation of electrocardiograph AsCARD B5 ECO Mr GREEN is made with the help of functional keyboard and LCD display which are placed in the upper part of the unit.

To start the electrocardiograph switch on power supply with the use of  key. If you do not press any key for 3 min the unit switches off automatically. To switch the unit on again press the  key.

### 3. 2. Arrangement of the keys.

-  power supply ON key
-  power supply OFF key
-  automatic examination start key
-  manual examination start key
-  stop paper key

copy key

change of leads key

sensitivity key

1 or 3-channel registration key




filters key


speed key

### 3. 3. Signalisation on LCD.

- **5, 10 lub 20mm/mV** signalises the chosen sensitivity,
- **2.5, 5, 25 lub 50mm/s** signalises the chosen speed,
-  signalises the chosen automatic examination,
-  signalises the chosen manual examination,
-  signalises the chosen copy of automatic examination,
- **35Hz, 50Hz, 35/50Hz** signalises filters on/off,
- **INOP** signalises the poor contact of the electrodes with patient body. It is impossible to perform automatic examination till INOP is signalised. Start of examination record with INOP displayed gives straight lines in all 12 leads. It doesn't matter which electrode is badly connected. Check the contact on each electrode.
-  signalises out of paper in paper box,
-  signalises discharging of internal accumulator,
- **I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6** signalises the chosen lead or group of leads.

### 3. 4. How to put ECG paper into the unit.

Lack of paper is signalised with  mark. To change the paper in the unit, remove the cover of the paper box and take out the finished roll. Put new roll inside the box and slip the paper through the driving shaft. The unit should take the paper automatically. To help the paper slipping, press  key. Then put the cover of the paper box on.

To prevent the linear head from pollution and the unit from blocking, it is necessary to use thermoactive and dustfree paper, which is always available at Aspel.

**Technical Specification of ECG paper:**

- |  |                           |
|--|---------------------------|
| - Kind of paper                        | thermoactive, free of wax |
| - Kind and colour of overprint         | millimetre net, orange    |
| - Width of paper                       | 58 mm ± 1 mm              |
| - Diameter of a roll of paper          | 43 mm ± 1 mm              |
| - External diameter of the paper shaft | 16 mm ± 0,5 mm            |

Note !

The producer of electrocardiographs do not take the responsibility for any defects caused by using the other kind of paper.

**3. 5. Arrangements of the electrodes.**

The electrocardiograph AsCARD B5 ECO Mr GREEN is equipped with a 10-electrode patient's cable.

In order to receive 12 standard leads of Einthoven, Goldberger, Wilson, the electrodes should be attached in the following way:

**Limb electrodes:**


- |            |             |
|------------|-------------|
| R - red    | - right arm |
| L - yellow | - left arm  |
| F - green  | - left leg  |
| N - black  | - right leg |

**Precardiac electrodes :**

C1; white-red            the fourth intercostal space on the right side of the sternum

C2; white-yellow  
C3; white-green  
C4; white-brown  
C5; white-black  
C6; white-violet

the fourth intercostal space on the left edge of the sternum  
half way between C2 and C4  
the fifth intercostal space in the left central clavicular line  
straight line from the C4 point perpendicularly to the left front  
auxiliary line in the intersection point with this line.  
at the same level as C5, but in left central auxiliary line.

### **How to clean the electrodes.**

To clean the electrodes use a duster moistened with water and soap.

### **The electrodes have to be disinfected after every examination.**

If the patient has been simultaneously connected to the several units, one should remember to estimate any possible risks that may have the result from leakage current summing for each of used units.

Pay attention !

During the electrodes connection, the conducting parts of electrodes and patient cable should not be in contact with each other or with other metallic parts, including ground connection.

The electrocardiogram shall appear in no longer than 10s after defibrillator impulse.

### **Protection against defibrillating impulse.**

The input circuit of AsCARD B5 ECO Mr GREEN is protected against defibrillating impulse.

## 4. EXAMINATION PERFORMANCE

### 4. 1. How to perform manual examination.

Switch on the unit with the use of the  key. The display shows the current set registration parameters. The unit supervises the electrodes contact with the patient's body and signalises poor contact with **INOP**.

When **INOP** is signalised, the ECG leads are printed as straight lines. It doesn't matter which electrode is badly connected. Check the contact on each electrode.

When done:

Choose the record speed

Choose the record sensitivity

Choose the lead to record

Choose the registration mode

Set the filters if required

Press the  key. The unit starts electrocardiogram registration - the chosen ECG leads, speed, sensitivity and filters are presented on the display.

During the electrocardiogram registration it is possible to change such parameters as:

- ECG leads
- record speed
- record sensitivity
- filters switching on or off
- record mode

Every parameter's change is automatically described on the display and registered on the electrocardiogram.

ECG registration lasts until the  key is pressed or ECG paper is finished.

Press the  key the registration is broken off. Press the  key to switch off the unit.

NOTE !

If the unit is battery operated and no keys are pressed or nothing is printed, it switches off automatically after 3 minutes.

To switch the unit on press the  key.

### 4. 2. How to perform automatic examination.

Switch on the unit with the use of the  key. The display shows the current

registration parameters. The unit supervises the electrodes contact with patients body and signalises poor contact with **INOP**,

When **INOP** is signalised, the ECG leads are printed as straight lines. It does not matter which electrode is badly contacted. If the INOP is displayed the auto examination will not be carried out. Check the electrodes contact before the examination start.

Choose the record speed

Choose the record sensitivity

Choose the lead to record

Choose the registration mode

Set the filters if required

Press the  key. The unit stores the ECG in the memory simultaneously from all 12 leads. After memorising 5 sec course, the unit starts the ECG printout (with parameters described before the examination start). The length of ECG leads is described in the unit configuration.

The examination can be interrupted by pressing the  key.

Press the  key to switch off the unit.

#### **4. 3. How to perform copy of automatic examination.**

Switch on the unit with the use of  key. The display shows the current registration parameters.

Choose the record speed

Choose the record sensitivity

Choose the registration mode

Set the filters if required

Press the  key. If there is any correct auto examination in the memory the unit starts to print the memorised electrocardiogram in the mode described before the print start. The duration of ECG leads is 5 sec.

Examination can be interrupted by pressing the  key.

Press the  key to switch off the unit.

## 5. HOW TO CHANGE REGISTRATION PARAMETERS

After switching on the unit the current chosen registration parameters are presented on LCD display. Information showing on the display are printed on the paper as description of registered ECG leads. Every parameters change, during the registration, causes the print of ECG status on the printout.

The change of sensitivity, leads or filters causes additional printout of 1mV pattern.

### 5. 1. Record speed.

It is possible to record in 5; 25; 50 mm/s. The chosen value is presented on the display. To change the speed use the  key. Every change of speed, during the registration, causes the print of ECG status on the printout.

### 5. 2. Record sensitivity.

It is possible to set 2,5; 5; 10; 20 mm/mV record sensitivity. The chosen value is presented on the display. To change the sensitivity use the  key. Every change of the sensitivity, during the registration, causes the printout of 1mV pattern and present status.

### 5. 3. ECG leads.

To change the leads use the  key.

The leads are changing from I to V6. Every change of leads, during the registration, causes the 1mV pattern printout with the leads description.

### 5. 4. Additional filter.

It is possible to set additional filters of ECG signal:

50 Hz	mains disturbances filter
35 Hz	muscular disturbances filter
35 Hz/50 Hz	muscular and the mains disturbances filter

Every change of the filters, during the registration, causes the printout of 1mV pattern and current status. To switch on/off filters use the  key.

## 5. 5. Registration mode.

It is possible to make ECG examination in below given modes:

- |            |                               |
|------------|-------------------------------|
| 1 channel  | printout of 1 lead            |
| 3 channels | sequenced printout of 3 leads |

Every change of the mode, during the registration, causes the printout of 1mV pattern and current status. To change the modes use the  key.

## 6. ACCUMULATOR CHARGING

Discharging of internal accumulator is signalled on the display by  $\text{+|}$ . To avoid improper functioning of the unit, the accumulator has to be charged.

Connect external power supplier of direct current, output voltage 12-24 V, power 3 W, polarity as on description plate, nearby its socket (negative polarity from the middle of plug) to the socket placed at the back side of the unit.

### NOTE !

It is allowed to perform ECG examination with the external power supplier plug in. It will not cause a danger to a patient and operator.

It is forbidden to use the external power supplier that does not fulfill the above given parameters.

Completely discharged accumulator has to be charged for 12-14 hours.

### ATTENTION !

There is a battery containing the lead in the apparatus. Don't throw the used up battery to a rubbish - bin! Subject it to recycle according to the obligatory regulations.

## 7. HOW TO MAINTAIN THE UNIT

To clean the cover use a duster moistened with water and soap.

**DO NOT USE SOLVENT AND CORROSIVE LIQUID !!!**

## **8. KEEPING CONDITIONS**

Electrocardiograph AsCARD B5 ECO Mr GREEN has to be kept in below given conditions:

- temperature: +5 + 40 °C,
- relative humidity: < 80% (no condensation),
- the air should not be polluted with corrosion ingredients

## **9. ELECTROCARDIOGRAPH AsCARD B5 ECO Mr GREEN EQUIPMENT**

1. Limb electrodes (EKK)\* 4 pcs
2. Chest electrodes (EPP)\* 6 pcs
3. ECG Patient cable - 1pc
4. Paper 58 mm - 1 roll
5. ECG gel - 250 ml
6. Operation manual - 1 pc

\* the electrodes (type EKK and EPP) and other accessories are available at Aspel.

## 10. TYPICAL PROBLEMS OF THE UNIT.

<b>Problem</b>	<b>Cause</b>	<b>How to solve</b>
Paper doesn't run	Out of paper paper was improperly inserted Accumulator was discharged	Put new roll of paper press key and check if paper wasn't blocked Switch off the unit and charge accumulator
Unit switches off automatically	Unit switches off automatically once any key wasn't pressed for 3 min Accumulator was discharged Thermal printer was overheated	Correct operation  Switch off the unit and charge accumulator Check the disturbances, switch filters on
Straight line are printed in all leads	Improper contact between patient body and electrode/s (INOP)	Correct contact of electrode/s, add more ECG gel, check patient cable connection
Printout in poor quality	Discharge of accumulator Incorrect ECG paper	Charge accumulator Change ECG paper
Auto examination can't be performed	Improper contact between patient body and electrode/s (INOP)	Correct contact of electrode/s, add more ECG gel, check patient cable connection
High level of disturbances	Poor contact between electrode/s and patient body and/or disturbances come from other switching equipment and/or patient is moving during examination	Correct contact of electrode/s, add more ECG gel, switch on digital filters, switch off other equipment, ask patient not to move during examination
Printout exceeds size of ECG paper	Too high sensitivity	Decrease sensitivity

**TABLE OF CONTENTS**

1.	APPLICATION	23
1.1	Description of the unit.	24
1.2	Picture of the unit .	24
1.3	Basic Technical Parameters.	24
1.4	Manufacturer.	25
1.5	Marking of the unit.	25
2.	WORKING CONDITIONS.	25
3.	FEATURES.	25
3.1	How to prepare AsCARD B5 ECO Mr GREEN to work.	26
3.2	Arrangement of the keys.	26
3.3	Signalisation on LCD	27
3.4	How to put ECG paper roll to the unit.	27
3.5	Arrangement of electrodes.	28
4.	EXAMINATION PERFORMANCE.	30
4.1	How to perform manual examination.	30
4.2	How to perform automatic examination.	30
4.3	How to perform copy of automatic examination.	31
5.	HOW TO CHANGE REGISTRATION PARAMETERS.	32
5.1	Record speed.	32
5.2	Record sensitivity.	32
5.3	ECG leads.	32
5.4	Additional filters.	32
5.5	Registration mode.	33
6.	ACCUMULATOR CHARGING.	33
7.	HOW TO MAINTAIN THE UNIT.	33
8.	KEEPING CONDITION.	34
9.	ELECTROCARDIOGRAPH AsCARD B5 ECO Mr GREEN EQUIPMENT.	34
10.	TYPICAL PROBLEMS WITH THE UNIT.	35