ASPEL s.a.

32-080 ZABIERZÓW, os. Sienkiewicza 33
tel. +48 12 285-22-22, fax +48 12 285-30-30

www.aspel.com.pl

Elektrokardiograf

AsCARD Mr. Silver

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ZABIERZÓW, luty 2005

WYDANIE II

Mr. Silver

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EDITION II

Mr. Silver

- ***Congratulations on your purchase of AsCARD Mr. Silver which is the effect of our many year research and experience gathered thanks to close contacts with our customers. You have chosen high quality, durability and efficiency - features characteristic for AsCARD Mr. Silver***
- ***ASPEL offers a wide assortment of accessories for ECG units, such as trolleys, ECG bags, cables, electrodes and recording paper.***
- ***Read this Instruction carefully as it contains guidelines on safe installation, operation and maintenance, as well as certain practical advice which will help to optimise the manner in which the unit is used.***
- ***Store this Manual for reference.***

APPLICATION

The electrocardiograph AsCARD Mr. Silver three, six and twelve channels unit enables you to make electrocardiogram in full 12 leads. The ECG unit maintenance is very simple and the modern keyboard makes the operation much easier. The record is made in manual or automatic mode with the possibility to make analysis and interpretation. The unit can be supplied from the net with 90-250V voltage or from internal accumulators built-in the unit. The accumulator capacity makes possible to print about 60 examination in auto mode. Low magnitude and small size are the additional features which enable easy transport. It is possible to carry out spirometry examinations and to send ECG examinations by the phone. Electrocardiograph is designed for carrying out ECG examinations at the patients home and in all health-care institutions.

1. The person operating the electrocardiograph ought to read the Operation Manual carefully before starting the unit operation.
2. The Operation Manual will help the user to operate and maintain electrocardiograph correctly.
3. The efficient functioning and fulfilling the tasks the user reaches after obeying the Operation Manual descriptions.
4. It is forbidden to use the unit with a damaged patient cable.
5. Every repair ought to be made by Authorised Services.
6. The manufacturer is not responsible for any damages resulting from not obeying the Operation Manual.
7. Electrocardiograph used together with a cardiostimulator or other electric stimulator do not cause any safety hazard the of patient or the user.
8. The efficiency of the unit and accessories shall be periodically controlled
9. If there are any irregularities in functioning of the unit, one shall contact Authorised Service, every time when they appear.
10. For the purpose of long-term archivisation it is requested to make copy of printed examination (xerocopy, scan). Thermoactive paper is too sensitive for the variable environmental conditions. It can have the result in illegibility of printout kept for a long period
11. Be extremely careful during defibrillation. It is prohibited to touch a defibrillated patient and equipment connected to the patient.


Electrocardiograph must not be used together with surgical devices working on high frequency.

1. DESCRIPTION OF ELECTROCARDIOGRAPH

1.1 Description of the unit.

The device is made up of modern microprocessor technology and all components are produced by the most reliable international manufacturers. The unit has a thermal printer, equipped with a head of high resolution to print the ECG leads. The cover made out of ABS resins what, together with the membrane keyboard facilitates cleaning.

1.2 View of the unit.

1. Paper box
2. RS 232C socket
3. Power supply socket
4. On/off switch
5. LCD display
6. Keyboard
7. Patient cable socket

1.3 Basic Technical Parameters.

Dimensions:	(LxWxH) 310x230x66 mm;
Weight:	2,3 kg;
Power supply:	-AC; 90 V÷250 V; 50/60 Hz ± 2%, -battery 12 V/1,3 Ah (built-in battery with the charger). Only authorised service centres are allowed to replace batteries; Battery capacity is sufficient for about 60 examinations;
Power consumption:	max 60 VA;
Working temperature: :	+10°C ÷ +40°C;
Relative humidity:	25% ÷ 95%;
ECG signals:	12 standard leads; I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6;
Sensitivity:	2,5/5/10/20 mm/mV ± 5%;
Recording speed:	5/10/25/50 mm/s ±5%;
Amplifier input impedance	>10 MΩ;
Paper:	thermoactive, dustfree, 112 mm wide;
CMRR:	> 100 dB;
Frequency range:	0,05 ÷150 Hz;
Sampling frequency:	1000 Hz;
Resolution:	2,8 µV/bit;
Skew between channels:	10 µs;
Digital filters:	50 Hz, 35 Hz i 25 Hz; anti-drift filter;
LCD:	graphical 320 x 240 pixels;
Usage safety:	- protection type CF (EN 60601-1) - protection class: I (EN 60601-1)
Class and group of the unit:	- class A, group 1 (CISPR-11)

ECG input circuit is protected against defibrillating impulse. ECG course will appear to 10 seconds after the defibrillating signal.

The electrocardiograph AsCARD Mr. Silver is the unit belonging to the A class. If used in households the unit may cause radio-waves disturbance. If the disturbance occurs the user may be forced to take precautions and bear responsibility for them.

Fuses:

- mains line **T800mA/L/250V**
- neutral line **T800mA/L/250V** (permitted to be replaced by the operator)
- internal power supply **B1 T6,3/L/250V** (to be replaced by authorized services only)

1.4 Manufacturer.

ASPEL S.A.
os. H. Sienkiewicza 33
32-080 Zabierzów / POLAND
tel. +4812 285 22 22 fax +4812 285 30 30

2. OPERATION, TRANSPORT AND MAINTENANCE CONDITIONS

Working conditions of electrocardiograph AsCARD:

- temperature: +10 to + 40 °C
- relative humidity: 25% to 95% (non-condensive)
- atmospheric pressure: 70 – 106 kPa

Transportation and storage conditions:

- temperature: -20 to + 60 °C
- relative humidity: 25% to 95% (non-condensive)
- atmospheric pressure: 70 – 106 kPa

Air should not be polluted with corrosioning redients.

3. INSTALATION OF ELECTROCARDIOGRAPH

After opening the packaging it is advised to check if equipment is complete. In case of any inconsistency please contact the seller or manufacturer.

After unpacking the device it is necessary to wait for approximately 1 hour, so the unit adjusts itself to climate conditions in the place it will be installed.

The unit has to be placed on a table, close to mains socket with 90...250V with frequency 50 or 60 Hz.

Connecting to the mains should be performed solely with original power supply cable delivered by manufacturer.

The device is not able to be used in places with combustibile gas presence or combustibile gas reeky.

It should be considered that in case of connecting the patient to several different devices the hazard resulted by total current has to be evaluated.

AsCARD is plugged in with cable to socket (3). The other end plug into AC 90...250 V_{AC} socket.

4. FUNCTIONS

The electrocardiograph AsCARD Mr. Silver allows the recording of 12 standard ECG leads. The following modes are possible:

- 3-channel manual record

this mode allows the recording of selected group of leads. It is possible to set required record speed, sensitivity, filters on/off. It is possible to print one of four groups leads I-II-III, aVR-aVL-aVF, V1-V2-V3, V4-V5-V6.

- 6-channel manual record

this mode allows the recording of selected group of leads. It is possible to set required record speed, sensitivity, filters on/off. It is possible to record one of two groups of leads - limb : I-II-III-aVR-aVL-aVF and precordial: V1-V2-V3-V4-V5-V6.

- 12-channel manual record

this mode allows to record 12 leads I-II-III-aVR-aVL-aVF-V1-V2-V3-V4-V5-V6. It is possible to set required: record speed, sensitivity, filters.

- automatic record

this mode allows the recording of 10 sec of ECG signal from all 12 leads in three, six channel in AsCARD Mr. Silver. It then allows an automatic analysis including time intervals, segment amplitude measurement and the calculation of the electric axes. The full report printout includes ECG signal record, averaged P-QRS-T complexes, measurements and calculations results. Record parameters are the same as in manual record.

- spirometry examination record

this mode allows (after connecting spirometry attachment SPIRO 31) the recording of the spirometry examination including Flow Volume Capacity.

- copy of automatic record

this mode allows the recording of copy of automatic examination stored in the unit memory.

It is possible to print a few dozen examples of the examination on 112 mm paper.

ECG is fit for running examinations directly on the heart.

4.1 How to prepare the unit for operation.


Connect the electrocardiograph AsCARD to the power source.

Put the main cable to the mains socket (3) placed on the back side the unit and to the power source 90-250 V. The unit is switched on/off with the use of the switch (1) placed on the back side of the unit. The internal accumulator starts to charge automatically once the unit has been connected to the mains and switched on with the use of on/off switch. To switch the unit on press  the key. The display shows:

AsCARD	Mr. Silver
ver. 1.xx	

The device is equipped with automatic power supply switch off while working on internal battery mode. The unit switches off after 3 min, if any key has been pressed.


View of the back side of the unit

1. The mains switch
2. Fuses
3. The mains socket
4. RS socket

4.2 How to put ECG paper into the unit.

To put the paper into the unit, remove the cover of the paper box (no. 1 on View of the unit) then press the unit power supply . OUT OF PAPER appears on the display. Put the roll inside the box and slip the paper through the driving shaft. The unit should take the paper automatically. Then put the cover of the paper box on.

To prevent the linear head from pollution and the unit from blocking, it is necessary to use thermo-sensitive and dustfree paper.

Note:

Paper recommended for use with your units is always available at ASPEL.

Technical Specification of ECG paper - R-A4:


- | | |
|---|---------------------------|
| - Type of paper: | thermoactive, free of wax |
| - Kind and colour of overprint: | millimetre net, orange |
| - Width of paper: | 112 mm +0/-2 mm |
| - Diameter of a roll of paper: | 45 mm \pm 1 mm |
| - External diameter of the paper shaft: | 16 mm \pm 0,5 mm |

The producer of electrocardiographs do not take the responsibility for any defects caused by using the other kind of paper.

-  Record sensitivity
-  Copy of last auto examination
-  Record speed
-  Selection of 3, 6 or 12-channel recording
-  Spirometry examination
-  Remote transmission
-  Menu - configuration mode
-  Muscle filter
-  Mains filter

Graphical LCD display presents all the information needed in order to operate the device. The drawing shows sample screen:

4.4 Arrangement of the electrodes.

The electrocardiograph AsCARD is equipped with a 10-electrode patient's cable. In order to receive 12 standard leads of Einthoven, Goldberger, Wilson, the electrodes should be attached in the following way:


Limb electrodes:

- R - red - right arm
- L - yellow - left arm
- F - green - left leg
- N - black - right leg

Precardiac electrodes :

- C1 white-red the fourth intercostal space on the right side of the sternum
- C2 white-yellow the fourth intercostal space on the left edge of the sternum
- C3 white-green half way between C2 and C4
- C4 white-brown the fifth intercostal space in the left central clavicular line
- C5 white-black straight line from the C4 point perpendicularly to the left front auxiliary line in the intersection point with this line.
- C6 white-violet at the same level as C5, but in left central auxiliary line.

Arrangement of precordial electrodes

During the electrodes connection, the conducting parts of electrodes and patient cable should not be in contact with each other or with other metallic parts, including ground connection.







5. EXAMINATION PERFORMANCE

The device enables to input of patient data such as:

surname	sex
first name	blood pressure
ID	name of person ordering the examination
age	comment no. 1
height	comment no. 2
weight	

In order to input patient data, operator has to press  key while being in STOP mode. Display shows:

SURNAME

Data should be inserted following messages shown on display. The key  enables deleting last character. Polish fonts are available when you press  and a letter. Figures are available while simultaneous pressing of  and figure. After pressing  key the inserted data are input. Then the device enables inserting the following parameter. Exit from data insertion mode is performed by pressing  key or any other functional one. Combinations of:  and  keys - deletes all current patient data,  and  keys - deletes current parameter.

5.1 How to perform manual examination.

Switch on the unit with the use of the  key. The display shows:

STOP	I - II - III
25 mm/s	10 mm/mV

The unit supervises the electrodes contact with the patient's body and signalizes poor contact with INOP.

INOP	I - II - III
25 mm/s	10 mm/mV


When INOP is signalized, the ECG leads are printed as straight lines. It doesn't matter which electrode is badly connected. Check the contact on each electrode. When done:

Choose the record speed .

Choose the record sensitivity .

Choose the lead to record .

Choose the registration mode .



Press the key . The unit starts electrocardiogram registration - the chosen ECG leads, speed, sensitivity and filters are presented on the display.

PRINT	I	II	III
25mm/s	10mm/mV		

During the electrocardiogram record it is possible to change such parameters as:

- ECG leads
- record speed
- record sensitivity
- filters switching on or off
- record mode

Every parameter's change is automatically described on the display and registered on the electrocardiogram.

ECG registration lasts until the  key is pressed. Press the  key - then registration is stopped. Press the  key to switch off the unit.

5.2 How to perform automatic examination.

Switch on the unit with the use of the  key.

The display shows:

STOP	I - II - III
25 mm/s	10 mm/mV

The unit supervises the electrodes contact with patients body and signals any poor contact with **INOP**.

INOP	I - II - III
25 mm/s	10 mm/mV

When **INOP** is signalled, the ECG leads are printed as straight lines. It does not matter which electrode is badly connected. If the **INOP** is displayed the auto examination will not be carried out. Check the electrodes contact before the examination start.

Choose the record speed

Choose the record sensitivity


Choose the lead to record

Choose the registration mode

Press the  key. The display shows:

AUTO	I - II - III
25 mm/s	10 mm/mV

The unit stores the ECG in the memory simultaneously from all 12 leads. After memorising 10 sec course, the unit starts the ECG printout (described before the examination start). The length of ECG leads is described in the unit configuration.

The examination can be interrupted by pressing the  key.

If in configuration parameters the auto examination with analysis was set or during data collecting the  key was pressed, the unit starts to calculate the ECG parameters. The display shows:

PLEASE WAIT
*** ECG ANALYSIS ***

After 20 sec. the unit starts to print the report including:

- P-QRS-T complex of the most numerous class with marking beginnings and ends of segments. The lower line shows the class size. The record speed is always 50mm/s.
- table with the values of segments' amplitude in all leads
- verbal interpretation of P-QRS-T complex's calculated parameters

The manufacturer made its best to achieve correct and accepted standards. The results of analysis can not however stand for the sole criteria of electrocardiogram description and physician evaluation is necessary.

After finishing the printout the unit goes back to the beginning position:

STOP I - II - III
25 mm/s 10 mm/mV

Press the  key to switch off the unit.

5.3 How to perform copy of automatic examination.


Switch on the unit with the use of  key. The display shows:

STOP I - II - III
25 mm/s 10 mm/mV

- Choose the record speed
- Choose the record sensitivity
- Choose the registration mode

Press the  key. If there is any correct auto examination in the memory, the display shows:

COPY I II III
25mm/s 10mm/mV

The unit starts to print the memorised electrocardiogram in the mode described before the print start. The duration of ECG leads is described in the unit configuration. Examination can be interrupted by pressing the  key.

If in configuration parameters the auto examination with analysis was set, then the unit starts to print its results. When the printout is finished the unit comes back to the start position. The display shows:

STOP I - II - III
25 mm/s 10 mm/mV

Press the  key to switch off the unit.

6. HOW TO CHANGE REGISTRATION PARAMETERS

When the unit is switched on, LCD displays currently selected recording parameters. Those are the values set before switching the unit off.

STOP I - II - III
25 mm/s 10 mm/mV

Information showing on the display is printed on the paper as description of registered ECG leads. Every parameters change, during the registration, causes the print of ECG status on the printout. The change of sensitivity, leads or filters causes additional printout of 1 mV pattern.

6.1 Record speed.

It is possible to record in 5; 10; 25; 50 mm/s. The chosen value is presented on the display. To change the speed use the  key.

6.2 Record sensitivity.

It is possible to set 2,5; 5; 10; 20 mm/mV record sensitivity. The chosen value is presented on the display. To change the sensitivity use the  key. Every change of the sensitivity, during the registration, causes the printout of 1mV pattern and present status.

6.3 ECG leads.

To change the leads use the  key. The leads are changing from I to V6. Every change of leads, during the registration, causes the 1mV pattern printout with the leads description.

6.4 Additional filter.

It is possible to set additional ECG signal:

50 Hz	mains interference filter
35 Hz	muscular interference filter

25 Hz	muscular interference filter
35/50 Hz	muscular and mains interference filter
25/50 Hz	muscular and mains interference filter

It is possible to choose 60 Hz instead of 50 Hz. It can be made in "Configuration mode" - "Mains filter" (Chapter 7).

To switch 25, 35 Hz filter press the  key, to switch 50 Hz press the  key.

6.5 Registration mode.

It is possible to make ECG examination in below given modes:

3-channel	printout of 3 leads
6-channel	printout of 6 leads
12-channel	printout of 12 leads

To change the modes use the  key. Currently selected registration mode may be described based on description of leads:

3-channel

STOP I - II - III
25 mm/s 10 mm/mV

6-channel

STOP I - II - III - aR - aL - aF
25 mm/s 10 mm/mV

12-channel

STOP 12 LEADS
25mm/s 10mm/mV

7. HOW TO CONFIGURATE THE UNIT PARAMETERS

Electrocardiograph AsCARD Mr. Silver has a battery supported memory where is defined: data describing the starting position of the unit (after switching the unit on), last automatic examination and report format.

The parameters relating to the configuration printout, isoline filtration and ECG examination type can be changed in configuration mode only.

The meaning of the keys are changed in configuration mode. Arrows placed on the

left side of the keys are active.

Options can be selected with  and  keys. Selection is confirmed by  key.

 key is used to cancel parameter's change and to leave configuration mode. To enter configuration mode use  key. LCD display shows:

CHOOSE LANGUAGE:
>OK ENGLISH

It is possible to operate the device using one of the language available in the unit configuration, for example: POLSKI, ENGLISH, FRANCAIS, HRVATSKI, CESKY, SLOVENSKY, TURKCE, etc.

To set the language use   keys. Choose the required one and set ">" mark in OK position, then accept it with the use of  key.

The display shows:

ENTER SURGERY NAME

Insert the "Surgery name" with the use of alphanumeric keyboard. The description is printed once automatic examination is performed.

After pressing  key, the display shows:

INOP CONTROL
>FULL SIMPLIFIED

In simplified INOP control the device checks if electrodes are connected in a correct way and sinalizes lack of connection with INOP message on the display. The filters which controls the isoelectric line position are not switched on. In this setting electrocardiogram appears within 15s since the last electrode connection control.

In full INOP control the unit automatically controls the isoelectric line position. If any of the electrodes is badly connected, the ECG recording is impossible. After correcting the connection, the electrocardiogram appears within 3s since the last electrode was controlled.


To set the option use   keys. Set ">" mark in required position and accept with

 key, then next parameter is shown.

AUTO MODE WITH ANALYSIS:
> YES NO

If **YES** is selected, in automatic mode, the performed printout of all ECG leads will be followed by automatic analysis of memorised leads and the calculated results will be printed out.

If **NO** is selected, in automatic examination mode, the ECG will be printed and the unit will return to the start position. To run analysis press .

To set the option use   keys. Set ">" mark in required position and accept with  key. Next parameter appears on the display:

ECG INTERPRETATION:
> YES NO


Option **YES** - after calculation and table printout, the interpretation is printed. **NO** option switches interpretation off. To set the option use   keys. Set ">" mark in required position and accept with  key. The unit goes to the next parameters.

AUTO MODE PRINTOUT
5 [s]

Insert the required value from the keyboard (setting range, printing time: 0 to 10 seconds and accept with  key. The display shows:

ISOLINE FILTRATION:
> YES NO

Option **YES** activates ECG leads filtration.
Option **NO** switches off the filtration.

To set the option use   keys. Set ">" mark in required position and accept with  key, then next parameter is shown:

MAINS FILTER
> 50 Hz 60 Hz

The choice defines the frequency of mains interference filter. Set the value relevant to power supply system frequency. To set the option use keys. Set ">" mark in required position and accept with key, then next parameter is shown.

PRINTOUT PRECISSION
> 1 2 3

To set the option use keys. Set ">" mark in required ECG line thickness and accept with key, then next parameter is shown.

TIME 09 : 45

To set the time use the keys on the keyboard. Set hour, press "space" then set minutes. To accept the choice press .
The display shows:



DATE : 98 : 12 : 11

The date is set in the same way as time. Press key. The display shows next parameter:

REMOTE UNIT
> OFF MODEM


Selecting **OFF** option switches off communication via modem.
Press to open a window for selection of printing format.

FORMAT:
> STANDARD CABRERA

Select the required printout format with  . Accept with .

Option **PC** allows selection of a program for communication with the unit via modem or a mobile phone. The display shows:

REMOTE SOFTWARE
> CARDIOTEKA CARDIOTEL

Selection of CARDIOTEKA followed by pressing  will connect the unit with CardioTEKA program and will automatically send the examinations to PC. CARDIOTEL option allows selection of transmission mode:

TRANSMISSION ONLINE:
> YES NO


If you choose YES option for "Transmission online", after pressing  key, the unit will continually transmit ECG.

If you choose NO option and press  key - the unit will send the copy of automatic examination.


To switch on/off the display highlighting press  key. The display shows:

HIGHLIGHTING:
> ON OFF

To choose the required option use  . Accept with .

Pressing  once again enables the regulation of display brightness.

BRIGHTNESS:
> OK 50%

Use   (from 0 to 100%) to adjust the brightness. Accept the choice with  key. The unit goes to initial position. The display shows:

STOP I - II - III
25 mm/s 10 mm/mV

8. HOW TO PERFORM SPIROMETRY EXAMINATION

Connect the spirometry attachment SPIRO 31 cable into the patient cable socket (see picture on p. 38). Switch on the unit with the use of  key. The display shows:

STOP I - II - III
25 mm/s 10 mm/mV

To obtain spirometry option press the  key. The display shows:

	SEX	AGE	HEIGHT
PATIENT:	>M	40	168



Before starting the examination enter patient's data:


surname	sex
first name	blood pressure
ID	doctor's name
age	* 1: comment
height	* 2: comment
weight	

Press  key to enter patient's data. The display shows:

SURNAME

Data should be inserted following messages shown on display. The key  enables deleting last character. Polish fonts are available when you press  and a letter. Figures are available while simultaneous pressing of  and figure. After pressing  key, the inserted data are input. Then the device enables inserting the



following parameter. Exit from data insertion mode is performed by pressing  key or any other functional one. Combinations of:  and  keys - deletes all current patient data,  and  keys - delete current parameter.

Press  key. The display shows the set parameters.

EXPIRE	SEX	AGE	HEIGHT
PATIENT:	>M	40	168

A patient breathes out and then the unit automatically prints an examination report containing:

- flow curve - volume,
- volume curve in time,
- FVC - Forced vital expiratory capacity
- FEV 1,0 - One second forced expiratory volume
- PEF - Peak expiratory flow
- MEF 25,
- MEF 50,
- MEF 75 - maximum expiratory flow at 25%, 50% and 75% FVC before the end of forced expiration.

If an examination is saved in the unit memory, press  (in spiro mode) to print a spiro examination copy. Press  again to return to ECG mode.

9. CLEANING, DISINFECTION AND MAINTENANCE

Before you start cleaning, disinfections and/or maintenance switch the unit off and unplug.

It is recommended to clean the cover, cables and electrodes with a soft cloth or sponge dampened with soap water or mild detergent.

Solvents or other aggressive agents must not be used.
In case the unit gets flooded contact the service providers for servicing the unit.

It is recommended to clean the unit at least once a week. In case of frequent operation, the unit has to be cleaned more often.

After an examination is finished, put electrodes into a container filled with a liquid removing albumin, e.g. Sekusept Pulver 2% + activator 0,5% for about 30 minutes for disinfection. Then take electrodes out of the container, rinse and dry.

Apert of periodical inspection of all the cables, sockets and electrodes, the unit does not require any maintenance.

10. ENVIRONMENTAL PROTECTION

Operation of ECG does not produce any wastes dangerous for the environment.

A worn and torn unit should be utilised. You may send it back to the manufacturer or to a firm utilising electronic equipment.

ECG units are equipped with „Pb” batteries. Such batteries cannot be thrown out to litter bins but must be utilised according to applicable regulations.

11. ACCESSORIES

Limb electrodes (EKK)*	4 pcs
Chest electrodes (EPP)*	6 pcs
ECG Patient cable (KEKG30)	1 pc
Main cable (SN-14 type)	1 pc
Paper 112 mm (R-A4)	1 pc
ECG gel	1 pc
Operation manual	1 pc
Guaranty card	1 pc

All accessories are always available at ASPEL Customer Service.

12. TROUBLESHOOTING

The unit does not switch on

No power supply, discharged internal battery.
Plug the unit into the mains socket (check the position of on/off switch).

The unit stops at printing

Discharged internal battery.
Plug the unit to charge the battery.

INOP is signalised

The record from mains supply is possible
One of the electrodes is not connected.
Check positions of all the electrodes.
If INOP is still signalised it means that the patient cable is damaged. Contact the service.

STROB ERROR is signalised

The unit was kept in too low temperature.
Wait until the temperature goes up to at least +5°C.

Out of paper

Insert a new roll of paper into the box and slip under the driving shaft.

LCD displays „|F”

Battery is discharged.
Charge the battery.

At request, the manufacturer will deliver circuit diagrams, lists of components, descriptions, calibration instructions and other information which will help the qualified technical staff when they have to repair parts defined as repairable.

13. MEANING OF SYMBOLS

Important manufacturer's notes

Note! See associated documentation


Date of production

Manufacturer's address

Defibrillation resistant application part of CF type

Fuse

Maximum number of identical packing which may be piled one on another, where "n" is the maximum

Transport packing must be protected against moisture

Indicates proper vertical position of a transport packing.

Packing content is fragile and must be carefully handled

Temperature range

TABLE OF CONTENTS

	APPLICATION	36
1.	DESCRIPTION OF ELECTROCARDIOGRAPH	38
1.1	The description of the unit.	38
1.2	View of the unit.	38
1.3	Basic Technical Parameters.	39
1.4	Manufacturer.	40
2.	OPERATION, TRANSPORT AND MAINTENANCE CONDITIONS	40
3.	INSTALLATION OF ELECTROCARDIOGRAPH	40
4.	FUNCTIONS	41
4.1	How to prepare the unit for operation.	42
4.2	How to put ECG paper into the unit.	43
4.3	Arrangement of the keys.	44
4.4	Arrangement of the electrodes.	46
5.	EXAMINATION PERFORMANCE.	47
5.1	How to perform manual examination.	48
5.2	How to perform automatic examination.	49
5.3	How to perform copy of automatic examination.	50
6.	HOW TO CHANGE REGISTRATION PARAMETERS	51
6.1	Record speed.	51
6.2	Record sensitivity.	51
6.3	ECG leads.	51
6.4	Additional filters.	51
6.5	Registration mode.	52
7.	HOW TO CONFIGURATE THE UNIT PARAMETERS	52
8.	HOW TO PERFORM SPIROMETRY EXAMINATION	57
9.	CLEANING, DISINFECTION AND MAINTENANCE	58
10.	ENVIRONMENTAL PROTECTION	59
11.	ACCESSORIES	59
12.	TROUBLESHOOTING	60
13.	MEANING OF SYMBOLS	61

Mr. Silver