ASPEL S.A.

32-080 ZABIERZÓW, os. Sienkiewicza 33
tel. +48 12 285-22-22, fax +48 12 285-30-30
www.aspel.com.pl

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Ambulatory Blood Pressure Recorder

HoICARD CR-07

Operation Manual

ZABIERZÓW, February, 2004

EDITION II

CONGRATULATIONS

- ***We are pleased that you have chosen the HolCARD CR-07 Ambulatory Blood Pressure Recorder. Choosing the products offered by ASPEL you have proved that you do not approve compromise and you want the highest quality;***
- ***ASPEL has the pleasure to present the HolCARD CR-07 which is a result of years of surveys and experiences gathered through close contact with our customers. You have chosen high durability and efficiency as these are the features of HolCARD CR-07;***
- ***Read this Manual carefully as it includes guidelines concerning the safe installation, use and maintenance and several practical instructions how to optimise this apparatus application;***
- ***Store this Manual for reference.***

- A person operating a blood pressure recorder shall read carefully the Operation Manual and the Guarantee Card before taking it into operation.
- Operation Manual will help the user in correct operation and maintenance of the recorder.
- Complying with instructions provided by this Manual ensures efficient functioning of the unit.
- All repairs should be performed by authorized service points.
- The manufacturer shall not be held liable for any damages resulted from failing to observe instructions contained in this manual.
- Periodical check of the unit and accessories is required. All irregularities in functioning should be reported to an authorized service.

1. THE USE.

HoICARD CR-07 blood pressure recorder is designed for operation in clinics, hospitals, ambulatories, and surgeries. The unit registers blood pressure employing the oscillometric method. It is supplied with four LR6 (AA) alkaline batteries or (AA) accumulator enabling one week operation of the recorder or 500 measurements. The unit is made of subassemblies manufactured by the worldwide leaders in the advanced microprocessor technology. It is an ergonomic, light and small size unit.

2. FUNCTIONS.

HoICARD CR-07 allows one week operation or 500 measurements. Parameters are set (thanks to the implemented RS-232 interface) by the program installed on PC or by the recorder keyboard. The examination is stored in EEPROM until deleted by the user. LCD allows direct reading of the last recording and signals weak battery. EVENT button allows recording of the pressure on the patient's request.

3. BASIC TECHNICAL-OPERATIONAL PARAMETERS.

Size:	(LxWxH) 90 x 36 x 93 mm,
Weight:	<200 g,
Power supply:	four LR6 (AA) batteries or (AA) accumulators capacity of new alkaline batteries is sufficient for a week operation or 500 measurements,
Working temperature:	+10° C ÷ +40° C,
Relative humidity:	20% - 85% without condensation,
Atm. pressure:	700 hPa - 1060 hPa
Measurement method:	oscillometric
Measurement range:	0-300 mmHg (against atm. pressure),
Maximal pressure measurement error in the cuff for all the range: (static conditions):	±3 mmHg,
Maximal pressure measurement error for the whole range: (dynamic conditions):	±5 mmHg,
Usage safety:	BF protection type (PN-EN 60601-1), PC work- ing with a recorder shall meet requirements of PN-EN 60950
Protection class:	internally powered unit,
Fuse:	internal power supply - PolySwitch typu miniSM DE 190,
Usable characteristics:	one week operation or 500 recordings, built-in timer, weak battery (accumulator) alarm unlimited memory time, LCD, EVENT button.
Storage and transport:	
Temperature:	+5° C ÷ +40° C,
Relative humidity:	<80% without condensation,
Air pressure:	700 hPa - 1060 hPa

Note !

The recorder is not waterproof. Avoid damping and soaking.

4. MANUFACTURER AND REPRESENTATIVE.


Ambulatory Blood Pressure Recorder HolCARD CR-07 is manufactured by:

ASPEL S.A.
os. H. Sienkiewicza 33
32-080 Zabierzów, POLAND
tel. +48 12 / 285 22 22, fax +48 12 / 285 30 30

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
Cursor directions.

8. MEASUREMENT HANDLING.

Insert new LR6 (AA) batteries or charged accumulators (recommended capacity >950 mAh). Put a cuff round a patient's arm (see Cuff section) and go to recording mode via HolCARD software or press ↓ button to select an option:

START
YES

Define if the recorder is to cancel results of previous measurements ↓ .

DELETE
YES

Parameters can be changed with ← (NO) → (YES) buttons and activate recording mode by pressing ↓.

NOTE: In case of emergency, when the cuff is pumped but the air does not escape for several minutes, take the pipe off the recorder (letting out the air or taking the cuff off) and report the event to the recorder operator. Remember to check whether operation of the recorder does not worsen blood circulation in the patient (i.e. observe the limbs).

9. HOW TO READ THE RESULTS OF THE EXAMINATION.

Replace batteries if necessary, connect RS-232 cable to the recorder and read the result through the program installed in the PC (see instructions for use HolCARD CR-06). When the battery symbol is highlighted and other symbols faded, the batteries (accumulators) are discharged – replace batteries and read the result.

NOTE: If you use HolCARD CR-05 program, recorder must be activate before you read the results (flashing colon between hours and minutes). In case the unit is idle, press the EVENT or ↓ buttons once.

10. EVENT BUTTON.

EVENT button is used for changing a mode of recorder operation. When the recording is started (a pump starts operation) and a patient is not ready yet, press the EVENT button and the recorder will stop and restart in about 3 minutes. The patient will have time to get ready.

If a patient feels sudden change in arterial blood pressure, he/she may induce an additional measuring process by pressing EVENT button.

11. BATTERIES.

To replace batteries, remove the cover from battery compartment at the bottom of the recorder. Take the used batteries out and insert four new LR6 (AA) batteries according to the code situated on the casing. You may also use accumulators with capacity at least 950mAh. Remember about the process of self-discharging and charge accumulators directly before the use. The used batteries and accumulators shall be disposed in special containers according to environmental requirements. Weak batteries are signalled with highlighting of a battery symbol (other symbols are also highlighted). The completely discharged batteries are signalled with highlighted battery symbol (other symbols are dark), the unit ends the measuring process.

12. CHANGE OF RECORDING PARAMETERS.

Recording parameters may be set with the program installed in PC (see instructions for use of HolCARD program) or from the recorder keyboard. Press ↓ button to enter the parameter recording mode and the following information will be displayed:

START
YES

Next pressing ← button will display the following:

START
NO

Pressing ↓ button will move to setting of measurement periods.

The recorder enables to set 6 measurement periods with measurement frequencies. During setting the hours one has to remember that the periods decrease in 24-hours cycle successively one after another.

First period starts at time given by the user and lasts till the beginning of period 2. Last period lasts from the end of previous period till the beginning of period 1. To change the number of measurement periods the value --:-- (which responds the time 00:00) has to be set in the last period. While setting the time one have to remember that the time can change between previous and next period.


Change the time of period beginning with → (time forward) and ← (time back) with 5 minute steps.

When the required time is set, press ↓ button and go to setting the time of next period.


After setting all of the periods, set measurement frequencies in particular periods.

Change measurement frequencies with → (forward) and ← (back) with 5 minute steps. Accept with ↓ button and set next active measurement periods.

Change time (hours) with buttons → (time forward) and ← (time backward). Then set the date or return to the waiting mode with ↑ button:

The format of date is: DD:MM:YY. Set the date with buttons ← (back) and → (forward). When the date is set pass on to setting a language by pressing ↓ button (or return to the waiting mode with ↑ button):

**LANGUAGE
POLISH**

Select a language with → ← buttons. Then press ↓ to start setting measurement parameters (or go back to the waiting mode with ↑ button):

**MANOMTR
YES**

The buttons: (→) switches on and (←) switches off the display of the value of pressure in the cuff. Having selected the required parameter press ↓ to pass on to setting measurement parameter (or go back to the waiting mode with ↑ button):

**MEASURE.
YES**

The buttons: (→) switches on and (←) switches off the display of the result of measurement. Having selected the required parameter press ↓ to pass on to setting the deflation parameter (or go back to the waiting mode with ↑ button):

**DEFL.
3**

This parameter defines the speed of air escape from the cuff. Change values with ← to reduce the value and → to increase the value with 1mmHg step. Accept with ↓ button and pass on to working mode that depends on patient age.

**NEONATAL
NO**

Button → allows to switch on and button ← enables to switch off the mode "Neonatal". After choosing the mode ("Neonatal" - Yes) the recorder limits the maximal pressure

in the cuff up to 150 mmHg.

Accept the choice with ↓ button and go to setting the maximal pressure in the cuff. The value of the pressure depends on patient's age and expected maximal systolic pressure. The setting describes maximal pumping value to which the cuff can be pumped. This value has to be higher than expected pressure of about 30 mmHg. Maximal value of this pressure is 150 mmHg for "Neonatal" and 300 mmHg for other.

MAX PRES
300

Change the parameters with buttons: ← (decreasing) → (increasing) with 5 mmHg step.

Pressing ↓ button pass to the examination (or go back to the waiting mode with ↑ button). This option is designed for calibrating the manometer.

TEST
mmHg 0

Pressing → button pump the cuff. Pressure value is displayed. Let out the air with ← button. When the required recording parameters are set, press ↓ button to start the measuring process. The following will be displayed:

MEASURE.
YES

When press ↓ button again, the recorder will ask whether to delete previous results or to continue recording (i.e. after replacing of the batteries we want to save the previous results and continue recording).

DELETE
YES

Change parameters with → ← buttons and start measurement with ↓ button. As the recording process starts with activation of the pump and pressure check, remember to put the cuff on as described in Cuff section) and place the unit in a bag. If MANOMETER is set at YES, recorder will display cuff pressure value. Results of

the last measurement (if MEASUREMENT is set at YES) will be displayed for about 3 seconds after the measurement is completed. .

Keyboard is non-active in measurement mode until batteries are removed.

13. MAINTENANCE.

Clean the casing with a soft cloth dampened with soapy water.

WARNING!

DO NOT USE SOLVENTS AND OTHER AGGRESSIVE AGENTS

14. RECORDER ACCESSORIES.

1. Cuff - adult, medium size - 1 pc
2. Bag for recorder with belt - 1 pc
3. RS-232 (KRSH-09) cable - 1 pc
4. LR6 (AA) alkaline batteries - 4 pcs
5. HolCARD CR-06 programme for operating the recorder - 1 pc
6. Operation manual - 1 pc
7. Guarantee card - 1 pc

15. TROUBLESHOOTING.

TROUBLE	CAUSE/HELP
Error displayed on LCD	Pressure in the cuff does not go up upon pumping or pressure remains in the cuff for too long time. <i>Check whether the cuff is correctly put on, cuff cable inserted into the unit socket, not twisted or squeezed. Make sure if the cuff or cable are not damaged.</i>
The unit switches off by itself (battery symbol highlighted, other symbols dark)	The unit switches off automatically by itself if batteries are not fit for operation (correct action). <i>To read the results replace batteries.</i>
No connection with PC	Error in transmission of RS232 signals. <i>Check whether RS232 cable is correctly connected to the recorder and to the PC, check whether COM port has been correctly selected in the transmission program, replace the cable with a new one, check batteries.</i> <i>The unit is idle, press ↓ and try to communicate again.</i>
EVENT button problem	<i>Problem with the EVENT button</i> <i>Tell the patient how to operate the EVENT button (section 10), check whether the button is not jammed.</i>

16. CUSTOMER SERVICE.

In order to receive the help in diagnosing the problem with HolCARD CR-07 you may contact service at:
ASPEL S.A. os. H. Sienkiewicza 33, 32-080 Zabierzów Poland, tel. +48 12 285 22 22,
fax +48 12 285 30 30

Table of contents

1.	The use.	21
2.	Functions.	21
3.	Basic technical-operational parameters.	22
4.	Manufacturer and representative.	22
5.	Location of regulating parameters.	23
6.	Marking.	24
7.	Preparation of the unit.	24
7.1	Cuff.	24
8.	Measurement handling.	25
9.	How to read the results of the examination.	25
10.	EVENT button.	26
11.	Batteries	26
12.	Change of recording parameters	26
13.	Maintenance.	30
14.	Recorder accessories.	30
15.	Troubleshooting.	31
16.	Customer service.	31